

MBH[®]

ANALYTICAL LTD

HOLLAND HOUSE • QUEENS ROAD • BARNET • EN5 4DJ • ENGLAND • TEL: +44 (0)20 8441 2024 • FAX: +44 (0)20 8449 0810
E-mail: info@mbh.co.uk Website http://www.mbh.co.uk

SETTING UP SAMPLE DATA SHEET

SUS Material Information

Sample Identity: **168X PBSUS 7 (batch B)**

Material type: LEAD BABBITT

Form: CAST DISC (~50mm dia x 20mm thick)

Produced by: MBH Analytical Ltd

Analysed by: MBH Analytical Ltd

Intended use: Routine spectrometer calibration checks


Composition (%)

Sb	7	Sn	0.1	Cu	0.05
As	0.7	Bi	0.05	Cd	0.005
Ag	0.01	In	0.002	Te	0.02
Ni	0.01	Zn	(0.02)	Au	0.005
Tl	0.01	Se	0.01	S	0.001

Notes

- 1) Use of this sample should be limited to the day-to-day checking of calibration stability for OE spectrometers.
- 2) The analytical data presented above are not certified as accurate, as this setting up sample is not intended to be used as a reference material for instrument calibration.
- 3) This material has been checked for homogeneity, and the data above are applicable to the first 15mm depth of the sample. The remaining 5mm should be discarded.
- 4) The material to which this document refers is supplied subject to our general conditions of sale.

MBH ANALYTICAL LIMITED



on 25th March 2013