

SETTING UP SAMPLE DATA SHEET

SUS Material Information

Sample Identity: **168X PBSUS PM1 (batch B)**

Material type: LEAD with PRECIOUS METALS

Form: CAST DISC (~48mm diameter)

Produced by: MBH Analytical Ltd

Supplied by: MBH Analytical Ltd

Intended use: Routine spectrometer calibration checks for fire assay

Analysis (ppm)

Pt	55	Pd	20	Au	35
Rh	12	Ru	0.1	Ir	2
Ag	40	Fe	1	Bi	100
Cu	5	Ni	3	Co	0.5
Te	1	As	2	Sb	1
Tl	10	S	2	Pb	(balance)

Notes

- 1) Use of this sample should be limited to the day-to-day checking of calibration stability for OE spectrometers.
- 2) The analytical data presented above are not certified as accurate, as this setting up sample is not intended to be used as a reference material for instrument calibration.
- 3) The data above are applicable to the whole sample.
- 4) The material to which this document refers is supplied subject to our general conditions of sale.

on 3rd January 2011

MBH ANALYTICAL LIMITED _____

C Eveleigh